

Specification

Safranine O (C.I. 50240) for microscopy

A0427

Physical Description:	Solid
Product Code:	A0427
Product Name:	Safranine O (C.I. 50240) for microscopy
Specifications:	Assay: min. 80 % λ_{\max} . (50 % EtOH): 530 - 534 nm E 1 %/1 cm, λ_{\max} : 1180 - 1480 Loss on drying: max. 15 %
WGK:	2
Storage:	RT
Molecular Formula:	$C_{20}H_{19}ClN_4$
M:	350.88 g/mol
CAS:	477-73-6
EINECS:	207-518-8
CS:	32041300

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode
 DRESDEFF508 • Finanzamt Darmstadt 007 225 00192 • Register court Darmstadt HRB Nr. 7340