

Specification

TES for buffer solutions

A1084

Physical Description:	Solid
Product Code:	A1084
Product Name:	TES for buffer solutions
Specifications:	Assay (titr.): min. 99 % pH (1 %; H ₂ O): 4.0 - 5.0 Residue on ignition: max. 0.1 % Heavy metals (as Pb): max. 0.0001 % Melting point: approx. 225°C Solubility (0.1 M; H ₂ O): clear, colorless Water: max. 1 % A (1 cm/0.1 M in H ₂ O) 260 nm: max. 0.05
WGK:	1
Storage:	RT
Molecular Formula:	C ₆ H ₁₅ NO ₆ S
M:	229.25 g/mol
CAS:	7365-44-8
EINECS:	230-906-3
CS:	29221900

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode
 DRESDEFF508 • Finanzamt Darmstadt 007 225 00192 • Register court Darmstadt HRB Nr. 7340