


Specification

**TAE buffer (50X)**

**A1691**

<b>Physical Description:</b>	Liquid
<b>Product Code:</b>	A1691
<b>Product Name:</b>	TAE buffer (50X)
<b>Specifications:</b>	<p>pH (20°C; H<sub>2</sub>O): 8.5 ± 0.2</p> <p><b>Composition:</b></p> <p>EDTA-Na<sub>2</sub> · 2H<sub>2</sub>O: 18.61 g/L (0.05 M)</p> <p>Acetic Acid glacial: 60.05 g/L (1 M)</p> <p>Tris: 242.28 g/L (2 M)</p>
<b>WGK:</b>	1
<b>Storage:</b>	RT
<b>CS:</b>	38221900
	

**AppliChem GmbH**