


Specification

Luminol

A2185


Physical Description:	Solid
Product Code:	A2185
Product Name:	Luminol
Specifications:	Assay (Perchl. Ac.): min. 95 % $\lambda_{\max. 1}$ (0.1 M NaOH): 344 -350 nm $\lambda_{\max. 2}$ (0.1 M NaOH): 297 - 303 nm Water (K.F.): max. 1 %
Hazard pictograms	
WGK:	nwg
Storage:	RT under argon
Signal Word:	Attention
GHS Symbols:	GHS07
H Phrases:	H315 H319 H335
P Phrases:	P261 P280 P305+P351+P338 P403+P233 P405 P501

AppliChem GmbH

Specification

Luminol

A2185

Molecular Formula:	$C_8H_7N_3O_2$
M:	177.17 g/mol
CAS:	521-31-3
EINECS:	208-309-4
CS:	29339980
	

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 007 225 00192 • Register court Darmstadt HRB Nr. 7340