

Specification

Peptone from Casein (pancreatic digest)

A2210

| | |
|------------------------------|--|
| Physical Description: | Solid |
| Product Code: | A2210 |
| Product Name: | Peptone from Casein (pancreatic digest) |
| Specifications: | <p>Enzyme activity: not detectable</p> <p>Total N: approx. 13 %</p> <p>pH (2 %; H₂O; 20°C): approx. 7.0</p> <p>Nitrogen (NH₂): 3.0 - 4.0 %</p> <p>Sulfated ash: max. 17 %</p> <p>NaCl: max. 6 %</p> <p>Water: max. 6 %</p> |
| WGK: | 1 |
| Storage: | RT |
| CAS: | 91079-40-2 |
| EINECS: | 293-428-4 |
| CS: | 35040090 |
| | |

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode
 DRESDEFF508 • Finanzamt Darmstadt 007 225 00192 • Register court Darmstadt HRB Nr. 7340