

Specification

Polyvinylpyrrolidone (K15) BioChemica

A2258

| | |
|------------------------------|---|
| Physical Description: | Solid |
| Product Code: | A2258 |
| Product Name: | Polyvinylpyrrolidone (K15) <i>BioChemica</i> |
| Specifications: | <p>Assay (from N): 11.5 - 12.8 %</p> <p>Solubility (10 %; H₂O): clear, colorless to yellowish</p> <p>pH (5 %; H₂O): 3.0 - 5.0</p> <p>Water (K.F.): max. 5 %</p> |
| WGK: | 1 |
| Storage: | RT |
| M: | approx. 10000 g/mol |
| CAS: | 9003-39-8 |
| CS: | 39059990 |
| | |

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 007 225 00192 • Register court Darmstadt HRB Nr. 7340