

Specification

Polyvinylpyrrolidone (K30) BioChemica

A2259

Physical Description:	Solid
Product Code:	A2259
Product Name:	Polyvinylpyrrolidone (K30) <i>BioChemica</i>
Specifications:	Assay (from N): 12.0 - 12.8 % Solubility (10 %; H ₂ O): clear, colorless to yellowish pH (5 %; H ₂ O): 3.0 - 5.0 Heavy metals (as Pb): max. 0.001 % Water (K.F.): max. 5 %
WGK:	1
Storage:	RT
M:	44000 - 54000 (g/mol)
CAS:	9003-39-8
CS:	39059990

AppliChem GmbH