

Specification

HEPES for cell culture

A3268

Physical Description:	Solid
Product Code:	A3268
Product Name:	HEPES for cell culture
Specifications:	Pyrogen test: passes test Assay (titr.): min. 99.5 % pH (1 %; H ₂ O): 4.7 - 6.0 Sulfated ash: max. 0.1 % Water: max. 0.5 % Chloride: max. 0.005 % Sulfate: max. 0.05 % As: max. 0.0001 % Fe: max. 0.0005 % Pb: max. 0.001 %
WGK:	1
Storage:	RT
Molecular Formula:	C ₈ H ₁₈ N ₂ O ₄ S
M:	238.31 g/mol
CAS:	7365-45-9
EINECS:	230-907-9
CS:	29335995

AppliChem GmbH