

Specification

L-Tyrosine for cell culture

A3401

Solubility:	0.4 g/L (H ₂ O)
Physical Description:	Solid
Product Code:	A3401
Product Name:	L-Tyrosine for cell culture
Headline Comment:	• All amino acids from AppliChem are of non-animal origin!
Specifications:	Pyrogen test: passes test Assay (titr.): min. 99 % α 20°C/D; 5 %, 1 N HCl, calc. on dried subst.: -12.3° - -11.0° Other amino acids: max. 0.5 % Ammonium: max. 0.1 % Chloride: max. 0.02 % Sulfate: max. 0.02 % As: max. 0.0001 % Fe: max. 0.001 % Pb: max. 0.001 %
WGK:	1
Storage:	RT
Molecular Formula:	C ₉ H ₁₁ NO ₃
M:	181.19 g/mol
CAS:	60-18-4
EINECS:	200-460-4
CS:	29225000

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode
 DRESDEFF508 • Finanzamt Darmstadt 007 225 00192 • Register court Darmstadt HRB Nr. 7340