

Specification

**Agar for molecular biology**

**A3477**

<b>Physical Description:</b>	Solid
<b>Product Code:</b>	A3477
<b>Product Name:</b>	Agar for molecular biology
<b>Specifications:</b>	DNases/RNases/Proteases: not detectable Ash: max. 4.5 % Gel point: 30 - 36°C Gel strength (1.5 %): 750 - 1100 g/cm <sup>2</sup> pH (1.5 %; H <sub>2</sub> O): 6.0 - 8.0 (60°C) Gelatin: not detectable Starch: not detectable Water: max. 12 %
<b>WGK:</b>	1
<b>Storage:</b>	RT
<b>CAS:</b>	9002-18-0
<b>EINECS:</b>	232-658-1
<b>CS:</b>	13023100

**AppliChem GmbH**

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • [info.de@itwreagents.com](mailto:info.de@itwreagents.com) • [www.itwreagents.com](http://www.itwreagents.com)  
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode  
 DRESDEFF508 • Finanzamt Darmstadt 007 225 00192 • Register court Darmstadt HRB Nr. 7340