



BIPEX200

Description of the container

Plastic container, HDPE (High Density Polyethylene), 220 L, conductive, outside black / inside natural.



PLASTIC TIGHT HEAD DRUM (in accordance with DIN ISO 20848-2)				
Material drum body:	Inside: Hostalen Histif – 5431Z, natural Outside: Hostalen GM-9350-C, black			
Colour body outside:	black			
Colour body inside:	natural			
Fabrication:	Extrusion blow moulded with two blow moulded bungs in the top, flat and bulged coat, top with blow moulded "TOP" ring, electrical conductive outer layer, interior layer of virgin made, interior electronically derivation by conductive connecting stripe, in the bottom one handling hole and one well (Sump) for totally drained			
UN - marking VPA 6	UN 1H1/Y1.6/200/*/D/ BAM11800-Schütz1 * year of production			
Resistance against compression characteristic number min.	25.000 N 600 N/mm			
Product drawing No.	3-43234			
Testing and description acc.	VPA 3, 4 and 8			

BUNG CLOSURE SYSTEM (in accordance with DIN ISO 20848-3) BCS 70 x 6 (fill side)

female

BCS 56 x 4

female

		Plug		
		Туре	closed	closed
Dimensions	(mm)	Thread:	male	male
Height	935 ±5	Material	PP	PP
Diameter	Ø 581 ±3	Colour	white	white
BCS 70x6 inside dia.	Ø 57.3	Gasket	O-ring	O-ring
BCS 56x4 inside dia.	Ø 52.9	Material	EBA	EBA
		Colour	blue	natural
Nominal capacity	220.0 ltr.	Cap seal	S70 loose	S56 loose
Brimful capacity min.	221.0 ltr.	Material	PP	PP
Weight without plugs	8.70 +0.2/-0.05 kg	Color	red	white
Weight with plugs	8.76 +0.2/-0.05 kg	Closing Torque (Nm)	hand screwed	15 +5/-2 N/m

Closure types

Bung thread